

Technical data sheet

Wireless accessories

Full discretion

Oticon Remote Control 3.0 offers an intuitive, straightforward and discreet way to operate hearing aids.

With the remote control, the hearing aid user can perform the most common everyday adjustments without drawing attention to the hearing aids:

- Volume up/down
- Select program
- Mute hearing aids

Full discretion is offered, since the remote control can operate the hearing aids from a distance of up to 1.8 metres without a direct line of sight between the remote control and the hearing aids.

Easy to operate

The remote control buttons are clearly visible and easy to operate.

Every time a button is pressed, a notification is heard in the hearing aids. At the same time the remote control confirms the action with a visual indication and a discreet click that can be heard and felt.

Remote Control 3.0



Remote Control 3.0 is a discreet and elegant solution for all users of Oticon hearing aids with 2.4 GHz radio technology.



Included

- Remote Control 3.0
- Batteries (2 x AAAA)
- Instructions for Use

| Technical specifications | |
|---------------------------------------|--|
| Model | RC-3A |
| Operating conditions | Temperature: 5 °C to 40 °C Humidity: 15% to 93 %, non-condensing Atmospheric pressure: 700 hPa to 1060 hPa |
| Storage and transportation conditions | Temperature: -25 to 70 °C Humidity: 15% to 93%, non-condensing |
| Mechanical dimensions | 75.3 x 35 x 12 mm; weight 32.7 g (incl. batteries) |
| Compatibility | Oticon hearing aids with 2.4 GHz radio technology |
| Battery type | 2 x AAAA |
| Battery life at normal use | 12 months |
| Typical transmission range | Up to 1.8 m |
| Status indicator | Green LED |
| Frequency | Licence free 2.4 GHz (ISM band) |
| Emission power | Max. 10 dBm EIRP |