


9 user-friendly modules. Clean, modern design.

## The New MedRx Studio Software

One Great Software for all the MedRx Products

Simple to use and easy to customize for your clinic




# The New MedRx Studio Software

## 9 user-friendly modules


Incorporate the New  
*MedRx Studio Software*  
into your daily practice to  
deliver quality professional  
hearing healthcare.

This clean, modern  
design is easy to use and  
functions as your complete  
Audiological suite.




### **DUAL CHANNEL AUDIOMETRY**

Perform air conduction, bone conduction, masking and speech testing using the mouse and/or keyboard with ease.




### **REAL EAR MEASUREMENTS**

IEC/ANSI standard REM functions are provided for: REUR, REUG, REIG, REAG, REAR, REOR, REOG and RECD.



### **HEARING INSTRUMENT TESTING**


The IEC/ANSI standard tests are available in an automatic sequence or for select measurements.



LSM

## LIVE SPEECH MAPPING


Live Speech Mapping (LSM) a quick and successful approach to hearing instrument fitting.



IMP

## TYMPANOMETRY


This function provides clear graphic and numeric displays of tympanogram and acoustic reflex data.



HLS  
MHA

## HLS/MHA


Easily simulate hearing loss and hearing instruments for patients.



VO

## VIDEO OTOSCOPE

Capture, view and store digital images. Retrieve results easily for follow-up visits.



VRA

## iVRA

A new form of visual reinforcement audiometry. A2D+ Audiometer, Stealth Audiometer & ARC compatible.



TI

## TINNOMETER


Revolutionary Tinnitus Assessment.

# M

## STUDIO

Software

*The MedRx Studio Software is easy to use and equips you with excellent counseling tools.*



**MedRx**  
Good Things Come in Small Packages

**MedRx International**  
MAICO Diagnostics GmbH  
Sickingenstr. 70-71  
10553 Berlin, Germany  
Phone: +49 30 70 71 46 50  
Fax: +49 30 70 71 46 99  
Email: medrx-sales@maico.biz  
www.medrx-diagnostics.com

Australian Distributor

**SONIC**  
EQUIPMENT

1800 639 263  
info@soniceq.com  
soniceq.com

**MedRx Minimum Computer Specs:**

Windows® PC Compatible Computer, Intel™ i5 Dual Core, 2.0 GHz Or Better. 4 GB RAM. 20 GB Free Hard Drive Space. Available 2.0 USB Port. Windows 7, 8 Or 10 Professional (32 Or 64-Bit).

MX-STUDIO-INTDS-1 Rev. 2

